







### Platinum-cured peroxide-free silicone hose

- No cross-sectional narrowing between screw connection and hose transition
- All wetted materials can be tracked
- Suitable for CIP and SIP

Product variants described in the data sheet may differ from the product presentation and description.

#### Can be combined with

	<b>Type BBS-4S</b> Hose fittings	▶
	<b>Type BBS-03</b> Fitting with welded connection according to DIN 11864-1 and BBS sterile contour	▶
	<b>Type BBS-05</b> Clamp connection with welded connection according to DIN 11864-3 and BBS-Sterile contour	▶
	<b>Type BBS-06</b> Sterile contour flange connection according to DIN 11864-2 with welded connection	▶

#### Type description

The platinum-cured silicone hoses are available with or without reinforcing fabric. We deliver the silicone hoses pre-assembled with the corresponding reusable hose screw connections BBS-4S or sold per meter without hose screw connection.

## Table of contents

<b>1. General technical data</b>	<b>3</b>
<b>2. Approvals and conformities</b>	<b>3</b>
2.1. General notes .....	3
2.2. Conformity .....	3
2.3. Standards .....	3
2.4. Foods and beverages/Hygiene .....	3
<b>3. Materials</b>	<b>4</b>
3.1. Bürkert resistApp .....	4
<b>4. Dimensions</b>	<b>4</b>
4.1. Hose Type HH, without reinforcement .....	4
4.2. Hose Type HI, glass silk-reinforced .....	5
4.3. Hose Type HJ, aramid-reinforced .....	6
4.4. Hose Type HK, double aramid-reinforced .....	7
4.5. Hose Type HL, stainless steel-reinforced .....	8
<b>5. Product installation</b>	<b>9</b>
5.1. Installation notes .....	9
Temperatures .....	9
Steam .....	9
Pressure .....	9
<b>6. Ordering information</b>	<b>9</b>
6.1. Bürkert eShop .....	9
6.2. Bürkert product filter .....	9
6.3. Ordering chart .....	9

## 1. General technical data

Product properties	
Dimensions	Further information can be found in chapter "4. Dimensions" on page 4.
Materials	
Wetted	Platinum-cured silicone Phtalate-free, latex-free, ADI-free
Performance data	
Operating pressure	Further information can be found in chapter "4. Dimensions" on page 4.
Medium data	
Operating medium	Gases or liquids
Medium temperature	Max. + 160 °C
Process/Port connection & communication	
Port connection	See <b>datasheet Type BBS- 4S</b> ▶
Approvals and conformities	
Certificate	REACH
Foods and beverages/Hygiene	
FDA	Further information can be found in chapter "2.4. Foods and beverages/Hygiene" on page 3.
USP	Further information can be found in chapter "2.4. Foods and beverages/Hygiene" on page 3.
EC Regulation 1935/2004	Further information can be found in chapter "2.4. Foods and beverages/Hygiene" on page 3.

## 2. Approvals and conformities

### 2.1. General notes

- The approvals and conformities listed below must be stated when making enquiries. This is the only way to ensure that the product complies with all required specifications.
- Not all available variants can be supplied with the below mentioned approvals or conformities.


### 2.2. Conformity

In accordance with the Declaration of Conformity, the product is compliant with the EU Directives.

### 2.3. Standards


The applied standards which are used to demonstrate compliance with the EU Directives are listed in the EU-Type Examination Certificate and/or the EU Declaration of Conformity.

### 2.4. Foods and beverages/Hygiene

Conformity	Description
FDA	<b>FDA – Code of Federal Regulations (valid for variable code PL02, P03)</b> All wetted materials are compliant with the Code of Federal Regulations published by the FDA (Food and Drug Administration, USA) according to the manufacturer's declaration.
USP	<b>United States Pharmacopeial Convention (USP) (valid for variable code PL04)</b> All wetted materials are biocompatible according to the manufacturer's declaration.
	<b>EC Regulation 1935/2004 of the European Parliament and of the Council (valid for variable code PL01, P02)</b> All wetted materials are compliant with EC Regulation 1935/2004/EC according to the manufacturer's declaration.

### 3. Materials

#### 3.1. Bürkert resistApp



#### Bürkert resistApp – Chemical resistance chart

You want to ensure the reliability and durability of the materials in your individual application case? Verify your combination of media and materials on our website or in our resistApp.

[Start chemical resistance check](#)

### 4. Dimensions

#### 4.1. Hose Type HH, without reinforcement

**Note:**

- Please note the information regarding steam sterilisation in chapter **"5. Product installation" on page 9.**
- A proper gap-free connection feature can only be guaranteed in combination with a BBS- 4S hose union.
- The available roll length has a tolerance of ± 10 %.
- The hoses can also be ordered pre-assembled in conjunction with BBS 4S hose unions.



Nominal diameter	Inner Ø	Wall thickness	Minimum bend radius	Max. operating pressure (at 20 °C)	Minimum burst pressure	Packaging unit	Functional part parameter	Article no.
	[mm]							
<b>DIN</b>								
DN 04	5	2.5	30	2	4	50	IA	730131
DN 06	6	3.0	50	2	4	50	IB	730137
DN 08	8	3.2	60	1.5	3	50	IC	730143
DN 10	10	4.75	70	1	3	25	ID	730148
DN 15	12.7	4.75	80	1	2	25	IE	730156
<b>ISO</b>								
DN 08	10	4.75	70	1	3	25	ID	730148
DN 10	12.7	4.75	80	1	2	25	IE	730156
<b>ASME</b>								
DN 1/4"	5	2.5	30	2	4	50	IA	730131
DN 3/8"	8	3.2	60	1.5	3	50	IC	730143
DN 1/2"	12.7	4.75	80	1	2	25	IE	730156

DTS 1000148925 EN Version: R Status: RL (released | freigegeben | valide) printed: 18.12.2024

#### 4.2. Hose Type HI, glass silk-reinforced

**Note:**

- Please note the information regarding steam sterilisation in chapter “5. Product installation” on page 9.
- A proper gap-free connection feature can only be guaranteed in combination with a BBS- 4S hose union.
- The available roll length has a tolerance of ± 10 %.
- The hoses can also be ordered pre-assembled in conjunction with BBS 4S hose unions.



Nominal diameter	Inner Ø	Wall thickness	Minimum bend radius	Max. operating pressure (at 20 °C)	Minimum burst pressure	Packaging unit	Functional part parameter	Article no.
	[mm]	[mm]	[mm]	[bar]	[bar]	[m]		
<b>DIN</b>								
DN 04	5	2	30	16	32	50	IA	730133
DN 06	6	3	50	16	28	50	IB	730139
DN 08	8	3.2	60	14	26	50	IC	730145
DN 10	10	4.75	70	12	24	50	ID	730152
DN 15	12.7	4.75	80	10	20	25	IE	730160
DN 20	19	5.8	150	8	18	25	IF	730168
DN 25	25.4	6.35	180	6	15	10	IG	730174
DN 32	32	7	220	4	12	10	IH	730180
<b>ISO</b>								
DN 08	10	4.75	70	12	24	50	ID	730152
DN 10	12.7	4.75	80	10	20	25	IE	730160
DN 15	19	5.8	150	8	18	25	IF	730168
DN 20	25.4	6.35	180	6	15	10	IG	730174
DN 25	32	7	220	4	12	10	IH	730180
<b>ASME</b>								
DN ¼"	5	2.5	30	16	32	50	IA	730133
DN ⅜"	8	3.2	60	14	26	50	IC	730145
DN ½"	12.7	4.75	80	10	20	25	IE	730160
DN ¾"	19	5.8	150	8	18	25	IF	730168
DN 1"	25.4	6.35	180	6	15	10	IG	730174





















DTS 1000148925 EN Version: R Status: RL (released | freigegeben | validé) printed: 18.12.2024

### 4.3. Hose Type HJ, aramid-reinforced

**Note:**

- Please note the information regarding steam sterilisation in chapter “5. Product installation” on page 9.
- A proper gap-free connection feature can only be guaranteed in combination with a BBS- 4S hose union.
- The available roll length has a tolerance of  $\pm 10\%$ .
- The hoses can also be ordered pre-assembled in conjunction with BBS 4S hose unions.




Nominal diameter	Inner Ø	Wall thickness	Minimum bend radius	Max. operating pressure (at 20 °C)	Minimum burst pressure	Packaging unit	Functional part parameter	Article no.
	[mm]							
<b>DIN</b>								
DN 10	10	4.75	70	16	30	50	ID	730150 
DN 15	12.7	4.75	80	16	30	25	IE	730158 
DN 20	19	5.8	150	12	24	25	IF	730164 
DN 25	25.4	6.35	180	10	20	10	IG	730172 
DN 32	32	7	220	8	18	10	IH	730178 
DN 40	38	7	260	5	12	10	II	730186 
DN 50	55	7	450	5	12	10	IK	730199 
<b>ISO</b>								
DN 08	10	4.75	70	16	36	50	ID	730150 
DN 10	12.7	4.75	80	16	30	25	IE	730158 
DN 15	19	5.8	150	12	24	25	IF	730164 
DN 20	25.4	6.35	180	10	20	10	IG	730172 
DN 25	32	7	220	8	18	10	IH	730178 
DN 32	38	7	260	5	12	10	II	730186 
DN 40	45	7	350	5	12	10	IJ	730193 
DN 50	55	7	450	5	12	10	IK	730199 
<b>ASME</b>								
DN 1/2"	12.7	4.75	80	16	30	25	IE	730158 
DN 3/4"	19	5.8	150	12	24	25	IF	730164 
DN 1"	25.4	6.35	180	10	20	10	IG	730172 
DN 1 1/2"	38	7	260	5	12	10	II	730186 
DN 2"	45	7	350	5	12	10	IJ	730193 

#### 4.4. Hose Type HK, double aramid-reinforced

**Note:**

- Please note the information regarding steam sterilisation in chapter “5. Product installation” on page 9.
- A proper gap-free connection feature can only be guaranteed in combination with a BBS- 4S hose union.
- The available roll length has a tolerance of ± 10 %.
- The hoses can also be ordered pre-assembled in conjunction with BBS 4S hose unions.



Nominal diameter	Inner Ø	Wall thickness	Minimum bend radius	Max. operating pressure (at 20 °C)	Minimum burst pressure	Packaging unit	Functional part parameter	Article no.
	[mm]	[mm]	[mm]	[bar]	[bar]	[m]		
<b>DIN</b>								
DN 15	12.7	4.75	90	25	48	25	IE	743018
DN 20	19	5.8	150	20	48	25	IF	743019
DN 25	25.4	6.35	180	15	32	10	IG	743020
DN 32	32	7	220	13	30	10	IH	743021
DN 40	38	7	260	10	20	10	II	743022
DN 50	55	7	450	8	20	10	IK	743024
<b>ISO</b>								
DN 10	12.7	4.75	90	25	48	25	IE	743018
DN 15	19	5.8	150	20	48	25	IF	743019
DN 20	25.4	6.35	180	15	32	10	IG	743020
DN 25	32	7	220	13	30	10	IH	743021
DN 32	38	7	260	10	20	10	II	743022
DN 40	45	7	350	10	20	10	IJ	743023
DN 50	55	7	450	8	20	10	IK	743024
<b>ASME</b>								
DN 1/2"	12.7	4.75	90	25	48	25	IE	743018
DN 3/4"	19	5.8	150	20	48	25	IF	743019
DN 1"	25.4	6.35	180	15	32	10	IG	743020
DN 1 1/2"	38	7	260	10	20	10	II	743022
DN 2"	45	7	350	10	20	10	IJ	743023

DTS 1000148925 EN Version: R Status: RL (released | freigegeben | validé) printed: 18.12.2024

#### 4.5. Hose Type HL, stainless steel-reinforced

**Note:**

- Please note the information regarding steam sterilisation in chapter “5. Product installation” on page 9.
- A proper gap-free connection feature can only be guaranteed in combination with a BBS- 4S hose union.
- The available roll length has a tolerance of ± 10 %.
- The hoses can also be ordered pre-assembled in conjunction with BBS 4S hose unions.



Nominal diameter	Inner Ø	Wall thickness	Minimum bend radius	Max. operating pressure (at 20 °C)	Minimum burst pressure	Packaging unit	Functional part parameter	Article no.
	[mm]							
<b>DIN</b>								
DN 10	10	4.75	70	16	35	25	ID	730154
DN 15	12.7	4.75	90	16	30	25	IE	730162
DN 20	19	5.8	150	12	24	25	IF	730170
DN 25	25.4	6.35	180	10	30	10	IG	743917
DN 32	32	7	220	8	18	10	IH	730182
DN 40	38	7	260	8	10	10	II	743852
<b>ISO</b>								
DN 08	10	4.75	70	16	45	25	ID	730154
DN 10	12.7	4.75	90	16	30	25	IE	730162
DN 15	19	5.8	150	12	24	25	IF	730170
DN 20	25.4	6.35	180	10	20	10	IG	743917
DN 25	32	7	220	8	18	10	IH	730182
DN 32	38	7	260	8	14	10	II	743852
DN 40	45	7	350	7	12	10	IJ	730195
<b>ASME</b>								
DN 1/2"	12.7	4.75	90	16	30	25	IE	730162
DN 3/4"	19	5.8	150	12	24	25	IF	730170
DN 1"	25.4	6.35	180	10	20	10	IG	743917
DN 1 1/2"	38	7	260	8,	14	10	II	743852
DN 2"	45	7	350	5	12	10	IJ	730195

DTS 1000148925 EN Version: R Status: RL (released | freigegeben | validé) printed: 18.12.2024



## 5. Product installation

### 5.1. Installation notes

#### Temperatures

The operating temperatures listed in this datasheet are valid only for the flow medium. The maximum values must be reduced for higher service life requirements. The limit values specified are temperatures that still guarantee sufficient flexibility of the hose. If the ambient temperature is not + 20 °C, consider this deviation when determining the maximum effective operating temperature. The specified pressure indications are valid also only at + 20 °C, unless specified differently. The correction values for higher pressures must be determined accordingly.

#### Steam

The BBS system silicone hoses may be steam sterilised at + 135 °C and a pressure of 3.5 bar. We recommend a maximum of 1.5 hours at + 135 °C. Between several sterilisation processes, at least one hour at room temperature should pass to achieve a hose stabilisation. Steam affects the mechanical and volumetric properties of the silicone elastomer. We recommend a visual hose inspection after 150 hours of steam exposure.

#### Pressure

If no specific information is provided, the correction factors may be used as guidelines. Be aware of the ambient temperature which has a significant influence on pressure resistance if the ambient temperature differs substantially from the pre-defined room temperature (+ 20 °C). BBS system hoses are not suitable for pulsating pressure applications. Likewise, continuous steam applications are excluded. Contrary to pipelines, silicone hoses are made of flexible, elastic material which accommodates the respective reinforcement material used to achieve a specific pressure resistance. Practically all such elastic material which are subject to different degrees of material fatigue and/or service life limitation. Therefore, hoses are always to be inspected resp. changed at regular intervals.

## 6. Ordering information

### 6.1. Bürkert eShop



#### Bürkert eShop – Easy ordering and quick delivery

You want to find your desired Bürkert product or spare part quickly and order directly? Our online shop is available for you 24/7. Sign up and enjoy all the benefits.

[Order online now](#)

### 6.2. Bürkert product filter



#### Bürkert product filter – Get quickly to the right product

You want to select products comfortably based on your technical requirements? Use the Bürkert product filter and find suitable articles for your application quickly and easily.

[Try out our product filter](#)

### 6.3. Ordering chart

#### Note:

Further information can be found in chapter [“4. Dimensions” on page 4](#).