Type 8200

Direct Welding

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© 2013 Bürkert SAS

Operating Instructions 1304/0_EU-ML_00564992_Original_FR

1. ABOUT THIS MANUAL	4	7.2. Installation onto the pipe	10
1.1. Symbols used.....	4	7.3. Inserting the probe	11
1.2. Definition of the word "holder"	4	7.3.1. Operating mode	11
2. INTENDED USE	5	7.3.2. The pressure evacuation holes feature	12
3. BASIC SAFETY INFORMATION	5	8. MAINTENANCE AND TROUBLESHOOTING	12
4. GENERAL INFORMATION	6	8.1. Safety instructions	12
4.1. International contacts.....	6	8.2. Maintenance	13
4.2. Warranty conditions.....	6	8.3. Cleaning	13
4.3. Information on the Internet.....	6	8.3.1. Cleaning in place (CIP).....	13
5. DESCRIPTION	7	8.3.2. Cleaning the interior of the direct welding holder	13
5.1. Area of application	7	8.3.3. Care of o-rings.....	14
5.2. Available version	7	9. SPARE PARTS AND ACCESSORIES	15
5.3. General description	7	10. PACKAGING AND TRANSPORT	15
6. TECHNICAL DATA	8	11. STORAGE	15
6.1. Conditions of use.....	8	12. DISPOSAL OF THE PRODUCT	16
6.2. Dimensions	8		
6.3. Materials	9		
7. INSTALLATION, COMMISSIONING	9		
7.1. Safety instructions.....	9		

1. ABOUT THIS MANUAL

This manual describes the entire life cycle of the direct welding holder. Please keep this manual in a safe place, accessible to all users and any new owners.

This manual contains important safety information.

Failure to comply with these instructions can lead to hazardous situations.

- This manual must be read and understood.

1.1. Symbols used

DANGER

Warns against an imminent danger.

- Failure to observe this warning can result in death or in serious injury.

WARNING

Warns against a potentially dangerous situation.

- Failure to observe this warning can result in serious injury or even death.

ATTENTION

Warns against a possible risk.

- Failure to observe this warning can result in substantial or minor injuries.

NOTE

Warns against material damage.

- Failure to observe this warning may result in damage to the fitting or system.

Indicates additional information, advice or important recommendations.

Refers to information contained in this manual or in other documents.

→ Indicates a procedure to be carried out.

1.2. Definition of the word "holder"

The word "holder" used within this manual refers to the direct welding holder type 8200.

2. INTENDED USE

Use of the direct welding holder that does not comply with the instructions could present risks to people, nearby installations and the environment.

- The 8200 direct welding holder is a holder welded onto a pipe and into which an analytical probe is inserted.
- Use the holder in compliance with the specifications and conditions of commissioning and use given in the contractual documents, in these operating instructions and in the operating instructions for the device which is inserted into it.
- Safe and trouble-free operation of the holder depends on its proper transport, storage and installation, as well as careful operation and maintenance.
- Only use the holder as intended.

→ Observe any existing restraints when the holder is exported.

3. BASIC SAFETY INFORMATION

This safety information does not take into account:

- any contingencies or occurrences that may arise during installation, use and maintenance of the holders.
- the local safety regulations for which the operating company is responsible including the staff in charge of installation and maintenance.

Danger due to high pressure in the installation.

Danger due to high temperatures of the fluid.

Danger due to the nature of the fluid.

Various dangerous situations

To avoid injury take care:

- Prevent any unintentional power supply switch-on.
- To ensure that installation and maintenance work are carried out by qualified, authorised personnel in possession of the appropriate tools.
- Use the holder only if in perfect working order and in compliance with the instructions provided in the instruction manual.
- To replace defective parts immediately.
- Observe the general technical rules when installing and using the holder.
- Not to fix, assemble or maintain the holder in explosive atmospheres.

Various dangerous situations

To avoid injury take care:

- not to use fluid that is incompatible with the materials from which the holder is made.
- not to use the holder in an environment incompatible with the materials from which it is made.
- not to subject the holder to mechanical loads (by placing objects on top of it or by using it as a step, for example).
- not to paint the internal or external threaded parts.

NOTE

The holder may be damaged by the fluid in contact with.

- Systematically check the chemical compatibility of the component materials of the holder and the fluids likely to come into contact with it (for example: alcohols, strong or concentrated acids, aldehydes, alkaline compounds, esters, aliphatic compounds, ketones, halogenated aromatics or hydrocarbons, oxidants and chlorinated agents).
- Contact with salt solutions and concentrated acids can result in corrosion.

4. GENERAL INFORMATION

4.1. International contacts

You may also contact your local Bürkert sales office.

The addresses of our international sales offices are available on the internet at: www.burkert.com

4.2. Warranty conditions

The condition governing the legal warranty is the conforming use of the 8200 in observance of the operating conditions specified in these operating instructions.

4.3. Information on the Internet

You can find these operating instructions and technical data sheets regarding the type 8200 at: www.burkert.com

5. DESCRIPTION

5.1. Area of application

The 8200 direct welding holder is a holder welded onto a pipe and into which an analytical probe is inserted. It is developed especially for applications in biotechnology and pharmaceuticals with great demands for cleanliness.

5.2. Available version

The direct welding holder exists in a single version with order code: **561728**.

The certificate 3.1 is supplied with all the holders.

5.3. General description


Fig. 1: Description of the holder

6. TECHNICAL DATA

6.1. Conditions of use

Fluid temperature	-30 to 140°C max. Temperature limits may depend on the inserted probe. Refer to the relevant operating instructions. If the temperature ranges given for the holder and the inserted probe are different, use the most restrictive range.
Pressure class	16 bar max. Pressure limits may depend on the inserted probe. Refer to the relevant operating instructions. If the pressure ranges given for the holder and the inserted probe are different, use the most restrictive range.

6.2. Dimensions


Fig. 2: Dimensions [mm] of the direct welding holder

6.3. Materials

Part	Material
Steel frame	Stainless steel
O-ring	EPDM
Seal pusher	Stainless Steel

7. INSTALLATION, COMMISSIONING

7.1. Safety instructions

DANGER

Risk of injury due to high pressure in the installation.

- Stop the circulation of fluid, cut-off the pressure and drain the pipe before loosening the process connections.

Risk of injury due to high fluid temperatures.

- Use safety gloves to handle the holder.
- Stop the circulation of fluid and drain the pipe before loosening the process connections.

Risk of injury due to the nature of the fluid.

- Respect the prevailing regulations on accident prevention and safety relating to the use of hazardous products.

WARNING

Risk of injury due to nonconforming installation.

- Fluidic installation can only be carried out by qualified and authorised personnel with the appropriate tools.
- Observe the installation instructions for the measuring device inserted into the fitting or the adapter.

Risk of injury due to an uncontrolled restart.

- Ensure that the restart of the installation is controlled after any interventions on it.

! WARNING

Danger due to nonconforming commissioning.

Nonconforming commissioning may lead to injuries and damage the holder and its surroundings.

- Before commissioning, make sure that the staff in charge have read and fully understood the contents of the manual.
- In particular, observe the safety recommendations and intended use.
- The installation must only be commissioned by suitably trained staff.

! WARNING

Gases or liquids can leak out undetected through the seals or the screws.

- Check regularly the holder for leaks.

! WARNING

Before initial operation of the holder, or after a long out-of use period:

- Ensure that the seal is tight.
- Ensure that all parts are in working order.

7.2. Installation onto the pipe

! WARNING

Risk of injury if the recommendations on installation of the measuring device inserted into the holder are not observed (see the relevant operating instructions).

- Take account of the recommendations on installation of the inserted probe.


Fig. 3: Dimensions [mm] for the installation of the holder onto the pipe

7.3. Inserting the probe

The holder is intended only for mounting probes with a length of 120 mm.

Make sure the welded area has cooled down before inserting the probe.


- Ensure there is no damage to the probe or the holder.
- Check that all o-rings are in place in their appropriate grooves, and are free of damage.
- Insert the seal pusher.
- Screw the probe into the PG 13.5 thread at a locking torque between 2 and 3 Nm.

7.3.1. Operating mode

When a probe is screwed into the direct welding holder, the seal pusher presses the o-ring into the space between the steel frame and the probe shaft. This seals the probe in a manner that almost completely eliminates crevices and their associated risk of contamination.

When the probe is loosened, the o-ring relaxes again. This lowers the friction between the o-ring and the other parts of the assembly, and the process is no longer sealed.

If this takes place while the process is still under pressure, the pressure evacuation holes feature (see chap. [7.3.2](#)), automatically comes into use.

Tip of a direct welding holder with loose probe: the o-ring is not under pressure.

Tip of a direct welding holder with tightened probe: the o-ring, under pressure, creates a seal with no crevices.

Fig. 4: Tip of a direct welding holder with loose or tightened probe

7.3.2. The pressure evacuation holes feature

The direct welding holder is equipped with a safety mechanism, to prevent danger to service personnel.

If the probe is unscrewed while there is still pressure in the process, the o-ring between the probe and the steel frame of the holder loses its seal. In this case, process medium enters the holder and exits through the pressure evacuation holes before the probe is loosened entirely (thereby informing the personnel that there is still pressure in the process).

In this case, tighten again the probe until the process is no longer under pressure. The probe can be then safely removed.

8. MAINTENANCE AND TROUBLESHOOTING

8.1. Safety instructions

DANGER

Risk of injury due to high pressure in the installation.

- Stop the circulation of fluid, cut-off the pressure and drain the pipe before loosening the process connections.

Risk of injury due to high fluid temperatures.

- Use safety gloves to handle the holder.
- Stop the circulation of fluid and drain the pipe before loosening the process connections.
- Keep all easily flammable fluid or material away from the holder.

Risk of injury due to the nature of the fluid.

- Respect the prevailing regulations on accident prevention and safety relating to the use of aggressive fluids.

WARNING

Risk of injury due to nonconforming maintenance.

- Maintenance must only be carried out by qualified and skilled staff with the appropriate tools.
- Ensure that the restart of the installation is controlled after any interventions.

8.2. Maintenance

WARNING

After repairs or maintenance.

- Ensure that the seal is tight.
- Ensure that all parts are in working order.

NOTE

Special attention must be given in the maintenance plan to the o-rings. The o-rings used in direct welding holder are consumable parts, and their life expectancy depends on the conditions in which they operate. If the temperature frequently exceeds 100°C, the o-rings must be replaced each time the probe is unscrewed, because they are easily deformed under conditions of high temperature. If replacement of o-rings according to the plan is not followed, leakage can occur.

- Check and replace o-rings as required; always replace at least once a year.
- Check that the probe is tightly screwed into the holder. (Weekly check recommended.)
- Perform CIP cleaning as required by the type of process.
- Clean the holder by means of the holes, if contamination with medium.
- Thoroughly clean the interior of the holder and the seal pusher if the probe breaks.

8.3. Cleaning

NOTE

The holder may be damaged by the cleaning product.

- Clean the holder with a cloth dampened with water or a detergent compatible with the materials the cleaning is made of.
- Never use scouring agent containing hard particles.

8.3.1. Cleaning in place (CIP)

If the installation into which the probe is built is being cleaned, the parts of the holder and the probe that are normally in contact with the medium (now in contact with the cleaning liquid) can be cleaned in place, together with the rest of the installation.

8.3.2. Cleaning the interior of the direct welding holder

If the space between probe and steel frame in the interior of the holder becomes contaminated, it can be cleaned easily by fitting cleaning adapters (not supplied by Bürkert) to the (non-threaded) holes in the wall of the steel frame, and injecting a suitable cleaning fluid. Alternatively, you can steam-sterilize.


Fig. 5: Cleaning of the holder with cleaning adapters



Fig. 6: Dimensions [mm] of a hole for cleaning adapters

8.3.3. Care of o-rings

- Check the o-rings from time to time.
- When mounting a new o-ring, take care not to damage the o-ring nut and the o-ring. Otherwise the tightness can no longer be guaranteed.

To replace the front o-ring:

- Forcefully pull the sleeve from the seal pusher
- Replace the o-ring

In higher temperature cycles (100 to 140°C) replace the o-ring each time the probe is unscrewed. This is because the compression set of EPDM must be considered (DIN ISO 815, the deviation of the original o-ring form over time, due to high temperature and pressure).

9. SPARE PARTS AND ACCESSORIES

ATTENTION

Risk of injury and/or damage caused by the use of unsuitable parts.

Incorrect accessories and unsuitable spare parts may cause injuries and damage the holder and the surrounding area.

- Use only original accessories and original spare parts from Bürkert.
- By default, the direct welding holder is delivered with EPDM o-rings.
- The adapter with order code **563477** is useful for mounting 120 mm oxygen or conductivity probes. The adapter is screwed onto the holder and the probe into the adapter. When in place, the end of the probe is almost flush with the top of the holder.


Fig. 7: Direct welding holder with adapter

10. PACKAGING AND TRANSPORT

NOTE

Damage due to transport

Transport may damage an insufficiently protected part.

- Transport the holder in shock-resistant packaging and away from humidity and dirt.
- Do not expose the holder to temperatures that may exceed the admissible storage temperature range.

11. STORAGE

NOTE

Poor storage can damage the holder.

- Protect the direct welding holder from cold, wet, direct heat sources (sun, heaters), soiling, and mechanical disturbance.
- Storage temperature : -30 to 140°C.

12. DISPOSAL OF THE PRODUCT

→ Dispose of the holder and its packaging in an environmentally-friendly way.

NOTE

Damage to the environment caused by products contaminated by fluids.

- Keep to the existing provisions on the subject of waste disposal and environmental protection.

Note

Comply with the national and/or local regulations which concern the area of waste disposal.

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