
## Type BBS-11

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## 1 THE OPERATING INSTRUCTIONS

The operating instructions contain important information.

- ▶ Read the instructions carefully and follow the safety instructions in particular.
- ▶ Keep the instructions in a location where they are available to every user.
- ▶ The liability and warranty for the product do not apply if the procedures in the instructions are not observed.

### 1.1 Symbols

→ designates a procedure which you must carry out.

**Warning of injuries:**

**DANGER!**

Immediate danger! Serious or fatal injuries.

**WARNING!**

Possible danger! Serious or fatal injuries.

**CAUTION!**

Danger! Moderate or minor injuries.

**Warning of damage:**

**NOTE!**

## 2 INTENDED USE

**Non-authorized use of the Type BBS-11 may be dangerous to people, nearby equipment and the environment.**

- ▶ Type BBS-11 has been designed as a safety ingold socket for sampling of equipment in the sterile area.
- ▶ Use according to the authorized data, operating conditions, and conditions of use specified in the contract documents and operating instructions.
- ▶ Correct transportation, storage and installation, as well as careful use and maintenance are essential for reliable and faultless operation.
- ▶ Use the product only as intended.

### 2.1 Restrictions

If exporting the products, observe any existing restrictions.

### 2.2 Definitions of terms

The term "product" used in these instructions always refers to Type BBS-11.

2

english

## 3 BASIC SAFETY INSTRUCTIONS

These safety instructions do not take into account any contingencies and events which may arise during the installation, operation and maintenance of the product.

**Danger – high pressure and discharge of medium!**

- ▶ Before loosening the nuts or screws, always turn off the pressure and relieve the lines/containers.
- ▶ Wear protective equipment if media is hazardous.

**General hazardous situations**

- ▶ Do not make any internal or external changes to the product.
- ▶ Ensure that the system cannot be activated unintentionally.
- ▶ Installation and maintenance work may be carried out only by authorized technicians with the appropriate tools.
- ▶ The product may be operated only when in perfect condition and in consideration of the operating instructions.
- ▶ As far as inspection, maintenance and repairs are concerned, observe national provisions of the country in which the connection elements are installed.
- ▶ The general rules of technology apply to application planning and operation of the product.

## 4 GENERAL INFORMATION

### 4.1 Contact address

**Germany**

Bürkert Fluid Control Systems  
Sales Center  
Christian-Bürkert-Str. 13-17  
D-74653 Ingelfingen  
Tel. + 49 (0) 7940 - 10 91 111  
Fax + 49 (0) 7940 - 10 91 448  
E-mail: info@de.buerkert.com

**International**

Contact addresses can be found on the final pages of the printed operating instructions. And also on the Internet at: [www.burkert.com](http://www.burkert.com)

### 4.2 Warranty

The warranty is only valid if the product is used as intended in accordance with the specified application conditions.

### 4.3 Information on the Internet

The operating instructions and data sheets for Type BBS-11 can be found on the Internet at: [www.burkert.com](http://www.burkert.com)

2

english

3

## 5 TECHNICAL DATA

### 5.1 Conformity

Type BBS-11 conforms to the EC directives according to the EC Declaration of Conformity.

### 5.2 Standards (if applicable)

The applied standards which are used to demonstrate compliance with the EC Directives are listed in the EC type test certificate and/or the EC Declaration of Conformity.

### 5.3 Identification

Information on material, pipe and connection dimensions can be found on the stamping on the product. The identification number of the product can be found on the supplied 3.1 certificate.


Fig. 1: Example of identification of the product

### 5.4 Sealing materials

| Seal material | Operating temperature                  |
|---------------|--|
| EPDM          | -40 °C to 90 °C, briefly up to 140 °C  |
| FEP           | -60 °C to 160 °C, briefly up to 205 °C |

Tab. 1: Sealing materials BBS-11

### 5.5 General technical data

|                                   |   |                                      |
|-----------------------------------|---|--------------------------------------|
| Material                          | comes into contact with medium                                    | Stainless steel 1.4435 BN2 (316L)    |
|                                   | does not come into contact with medium                            | Stainless steel 1.4401 or equivalent |
|                                   | Diaphragm   | silicone                             |
|                                   | O-ring  | EPDM/FEP                             |
| Pipe dimensions                   | see identification on product ("Fig. 1")                          |                                      |
| Permitted application temperature | Depending on sealing material, see "Tab. 1"                       |                                      |
| Ambient temperature               | -20 °C to +80 °C  |                                      |
| Media                             | Fluids  |                                      |
| Operating pressure                | -1 to +16 bar (depending on temperature and size, see data sheet) |                                      |

4

english

## 6 ASSEMBLY

### WARNING!

#### Danger – high pressure and discharge of medium.

- ▶ When working on the product or the system, always switch off the pressure and relieve the lines/containers.
- ▶ Wear protective equipment if media is hazardous.

#### Risk of injury from improper assembly.

- ▶ Installation must only be carried out by authorized technicians and with the appropriate tools.
- ▶ Secure system from unintentional activation.

### 6.1 Welding in the safety ingold socket

#### NOTE!

#### Leak due to damaged sealing elements.

- ▶ Do not weld in the assembled product. It is **essential** to remove the plug before welding and protect it from dust, flying sparks and other influences.

#### Leak due to damaged sealing contour.

- ▶ To ensure the sealing function, protect the sealing contour during installation, welding and cleaning procedures.
- ▶ When reworking the welded-in ingold socket, the sealing contour must not be damaged.

- ▶ The inner diameter  $D_i$  of ingold socket has to be machined to 25H7 respectively 40H7 after welding.
- ▶ Welding in and reworking of the safety ingold socket may be carried out only by adequately trained professionals.

→ Connect the parts **positively** in a protective gas shield.

When cleaning the weld seam by grinding or acid cleaning, observe the following before assembling the connection:

- Carefully remove all grinding dust and acid-cleaning residue.
- Do not damage the label.
- There must be no material abrasion on the sealing edges. Material abrasion will result in sharp-edged sealing contours and a damaged seal.
- Check sealing contour for damage.

We recommend preparing a welding report.

english

5

### 6.2 Assembly of dummy plug and Asepto sampling socket


Fig. 2: Assembly of dummy plug and Asepto sampling socket Type BBS-11

#### Dummy plug

- Ensure the correct position of the O-ring (2) on the dummy ingold (3). Screw the dummy ingold into the welded-in ingold socket (1) to lock it, using the union nut (4). Observe tightening torque (see "Tab. 2").

This way, sampling is no longer possible.

The retaining ring (5) is already installed upon delivery and must not be removed.

#### Asepto sampling socket

- Before screwing the ingold socket in, ensure the presence and correct position of the O-rings (2, 8, 10).

- Push the diaphragm (7) forward to the limit stop in the sampling socket (9) and screw the socket tightly to the sampling ingold (6). **Attention: This screw connection may not be loosened during operation or when under pressure.**

- Push this assembly into the welded-in socket (1) and screw down securely, using the union nut (4). Observe tightening torque (see "Tab. 2"). **Attention: This screw connection may not be loosened during operation or when under pressure.**

- Screw the locking element (11) onto the sampling socket. This is the only screw connection that may be opened during operation (for the purpose of sampling). The chain (12) prevents that the locking element gets lost.

6

english

|                        |            |    |
|------------------------|------------|----|
| Orifice (DN)           | 25         | 40 |
| Tightening torque (Nm) | hand-tight |    |

Tab. 2: Tightening torque of union nut Type BBS-11 ingold socket

### WARNING!

#### Danger – high pressure and discharge of medium.

- ▶ During operation and/or existing pressure, only the sampling cap (see "Fig. 2", sketch on the far right with chain) may be opened. It can be opened with just a few turns.
- ▶ All other screw connections may not be loosened in this state.

### 6.3 Disassembly

Disassembly is in reverse sequence to assembly.

## 7 START-UP

### WARNING!

#### Risk of injury from improper operation.

Improper operation may result in injuries as well as damage to the product and the surrounding area.

- ▶ Before start-up, ensure that the operating personnel are familiar with and completely understand the contents of the operating instructions.
- ▶ Observe the safety instructions and intended use.
- ▶ Only adequately trained personnel may start up the equipment/the product.
- ▶ Following assembly, ensure a controlled restart.
- ▶ When starting up the equipment, ensure that no unauthorized voltage increases and pressure surges can occur.

#### NOTE!

#### Replace wearing parts in time

- ▶ It is the responsibility of the plant operator to ensure the durability of the wearing parts. Observe that the wearing parts are replaced in time.

english

7

**NOTE!****Damage to the sealing elements when cleaning the system.**

- ▶ Clean the system preferably with a cleaning agent which does not damage the sealing elements.
- ▶ Do not clean with wire brushes or machines which cause abrasion of the surface.
- ▶ When using mechanical pipeline monitoring devices, ensure that they do not damage the sealing elements (also the sealing contour).

**Damaged sealing elements must be replaced.**

**8 MAINTENANCE****WARNING!****Danger – high pressure and discharge of medium!**

- ▶ Before loosening the nuts or screws, always turn off the pressure and relieve the lines/containers.
- ▶ Wear protective equipment if media is hazardous.
- ▶ Nuts or screws may be retightened on pressurized lines/containers only by technicians in consideration of special precautions.
- ▶ When shutting down the equipment, ensure that no unauthorized voltage increases and pressure surges can occur.
- ▶ Following maintenance, ensure a controlled restart.

**Have the product serviced regularly by technicians!**

Inspection and maintenance work includes in particular monitoring and ensuring the

- leak-tightness,
- identification,
- proper mode of operation of the safety and warning devices.

8

english

**9 SPARE PARTS****CAUTION!****Risk of injury and/or damage by the use of incorrect parts.**

Incorrect accessories and unsuitable spare parts may cause injuries and damage the product and the surrounding area.

- ▶ Use only original accessories and original spare parts from Bürkert.

Spare parts are: O-rings, silicone diaphragm.

The spare parts can be ordered from your Bürkert sales office. The order numbers for the spare part sets can be found on the corresponding data sheet on the Internet.

**NOTE!****Replace wearing parts in time.**

- ▶ It is the responsibility of the plant operator to ensure the durability of the wearing parts. Observe that the wearing parts are replaced in time.

**10 TRANSPORTATION, STORAGE, DISPOSAL****NOTE!****Transport damage.**

Inadequately protected sealing elements may be damaged during transportation.

- ▶ Transport the product in a firmly assembled state, protected against moisture and dirt, in shock-resistant packaging.

**Incorrect storage may damage the product.**

- ▶ Prevent the temperature from exceeding or dropping below the permitted storage temperature.
- ▶ Store the product in a dry and dust-free location!
- ▶ Store O-rings and diaphragm dry and protected from UV radiation and for not longer than 3 years.
- ▶ Storage temperature -40 to +80°C.

**Damage to the environment caused by products contaminated with media.**

- ▶ Dispose of the product and packaging in an environmentally friendly manner.
- ▶ Observe applicable disposal and environmental regulations.

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© Bürkert Werke GmbH & Co. KG, 2011 - 2017